

PEDIATRIC POWER BASE

The Quickie Z-500 provides the optimum power wheelchair base for pediatric seating and positioning making it ideal for the young rider. The unique seat and base design allows for multiple growth adjustments. Seat growth kits are also available. This model allows young riders to explore environments independently.

Z-500



Option



QUICKIE Z-500 & KID KART® XPRESS SEATING
Can be configured to accept Kid Kart Xpress seating for the easy switch between manual and powered mobility.



QUICKIE Z-500 POWER TILT
Get up to 30° of power tilt with our TK Power Tilt System with through drive control option.

